

permobil
POWER WHEELCHAIRS

READY OR NOT...



Koala
MINIFLEX




Dynamite performance in a small package means you can safely go wherever the adventure takes you.



MINIFLEX POWER SEATING


Unique power seat elevation and tilt designed to provide valuable positioning and exceptional low seat to floor height.



PARENTAL APPROVED

Optional features such as front bumper and remote stop give any parent that extra piece of mind for their little explorer.

COLOR OPTIONS



SPECIFICATIONS

Maximum User Weight	80 lbs
Maximum Speed	4 mph
Range	15 mi*
Base Width	23"
Base Length	32"
Minimum Turning Radius	21"
Ground Clearance	2.5"
Recommended Battery Type	M40
Weight Including Batteries	226 lbs
Drive Electronics	R-Net
Maximum Slope	10° Up, 6° Down†

Seat-to-Floor Height	15"
Optional Seat Elevator	8"
Power Tilt	0° – 25°
Manual Recline	85° – 120°
Manual Legrest Elevation	90° – 150°
Backrest Height	14", 16", 20"
Seat Depth	10", 12", 14"
Seat Width	10", 12", 14"
Armrest Pad Lengths	8", 10"

*Actual driving range will vary based on driving and battery conditions.
 †Refer to Owner Manual for safe driving practices.