

Over curbs, on trails, through grass and sand



You no longer need to spend time looking at the ground, worrying about every bump, stick or rock, or wasting energy wheeling over grass. With the FreeWheel you can look at the world around you and experience life to the fullest.

FreeWheel RACK Carry-All



With the RACK Carry-All you instantly become more independent and self-sufficient.

Weighing just 0.5 kg the RACK Carry-All works with your FreeWheel, allowing you to carry groceries or luggage when you go to the airport.

The RACK Carry-All measures 15" (38cm) x 20" (51cm) and easily "plugs in" to sockets that are attached to your wheelchair frame.

Will it fit my chair?

The FreeWheel is designed to work with most rigid-framed manual chairs. Our step-by-step Setup Guide shows how to install the FreeWheel for your specific footrest type and have you out and FreeWheelin' in no time.

To view the set up video, visit permobil.com.au.

Folding Wheelchair, no problem

The ADAPTOR for Folding Wheelchairs securely clamps on to the horizontal side tubing of your chair. Once set up for the individual chair, it comes on and off in seconds, making it easy to use. A standard FreeWheel attaches to the "hitch" out front instead of a wheelchairs footrest.



/PermobilAus



www.permobil.com.au



1300 845 483



info.au@permobil.com

available at

permobil

Freedom • Independence • Mobility

Permobil Australia © 2016 Permobil



FreeWheel

MOBILITY

FREEDOM

INDEPENDENCE



FreeWheel

www.permobil.com.au

What is the FreeWheel?

The FreeWheel is a lightweight, durable third-wheel that instantly turns your wheelchair into an all-terrain chair. This innovative adaptive equipment securely attaches to your solid footrest and safely lifts the front casters off the ground.

The FreeWheel was developed by Patrick Dougherty, following a motocross accident. While his wheelchair was crucial to his mobility, he quickly discovered it was also an obstacle to getting out and playing with his children.

As an engineer he knew the goal was to get the front casters off the ground. The end result, after years of testing, is the FreeWheel Wheelchair Attachment.

Now distributed around the world the FreeWheel is opening doors to wheelchair users that were once thought closed for good.

Whether it's participating in a race, at the park watching your child play sports, or camping and hitting the trails, the FreeWheel is the one piece of adaptive equipment you don't want to be without.



Super easy to put on, the FreeWheel transforms my chair into an all terrain chair in seconds so I can go where I want to. The FreeWheel is awesome!

Mal, Australia



My goal is to use my FreeWheel to explore all the wonders of the world!

Catherine, Canada



Rough terrain, mounds and bumps, Freewheel is now the reason I live life on the edge.

Paul, United Kingdom